ES USER'S GUIDE FOR THE PRACTITIONER

COMPOSITION

One glass vial containing: cellulose derivative product dissolved in ethyl alcohol to form a gel to which a contrast agent (tungsten) has been added.

PRESENTATION

DISCOGEL® comes in the form of a kit ready for the use with a descriptive notice and the following characteristics:

One type I glass vial containing 2.2 ml of solution for injection composed of a mixture ethyl alco- hol, cellulose derivative product and tungsten	Gelscom		2797		
Two 1ml syringe	Merit Medical		2797		
One high-flow 19G needle	Becton Dickinson		0318		
Two 18G spinal injection needles	Becton Dickinson		0050		

INDICATIONS

Treatment of herniated discs at any level of spine.

PROPERTY

The 96% pure ethyl alcohol produces a local necrosis of the *nucleus pulposus*. Its action is mechanical via a dehydration of the turgescnt and protruding disc, which is compressing the peripheral nerves of the rachis, and causing extreme pain.

DISCOGEL® is a sterile viscous solution containing ethyl alcohol, cellulose derivative product and tungsten: radio-opaque element. The injection into the disc of DISCOGEL® treats herniated discs and relieves lumbago, radicular or lombo-radicular. DISCOGEL® is available in the form of a kit containing 2.2 ml of solution for injection with two disposable 1ml syringes (each for single use) and 2 needles which make administration straightforward, assure sterility and prevent contamination during use. The quantity of gel in one kit (a single vial) permits filling the two 1 ml syringes provided.

INSTRUCTIONS AND DOSAGE

The medical act corresponds to an injection into the disc under radiological control. Injection is accomplished strictly into the herniated intervertebral disc whether it is lumbar, cervical or thoracic.

The resuspension of the metallic contrast agent must be homogeneous at the time of its injection❶. Kit allows intervention on two discs (one needle and one syringe per disc)❷. The quantity of DISCOGEL® injected per disc varies between 0.2 and 0.8ml, according to the dimension of the disc and extent of the hernia. A syringe of 1ml is provided so that any of these injection volumes is possible.

In general, it is recommended to use:

- 0.2 ml of DISCOGEL® for cervical discs
- 0.3 - 0.5 ml of DISCOGEL® for thoracic discs
- 0.6 - 0.8 ml of DISCOGEL® for lumbar and lumbar sacrum discs

PRECAUTIONS TO BE TAKEN DURING ADMINISTRATION

Injection

Act takes place in an operating room equipped for radiology and surgery under aseptic conditions. It is recommended that DISCOGEL® be injected at room temperature (see chapter "Precautions to observe"). DISCOGEL® is directly injected into the spinal disc.

The treatment is ambulatory.

Intervention

The injection into the disc is accomplished, after asepsis and local anesthesia via the post lateral❶ route for lumbar and thoracic discs and via the antero-lateral❷ route for cervical discs. Preferably the disc is punctured using a needle of:

- 20 Gauge for cervical disc (not provided in kit)
- 18 Gauge or thoracic and lumbar discs (provided in kit), so as to reach the central region of the intra spinal space

Injection of DISCOGEL® inside the disc

In absence of contraindication for injection of ethanol gel inside the disc, an injection intra-herniated disc of DISCOGEL® is accomplished into the nucleus pulposus. At the beginning of the injection the patient may experience a transiional scalding sensation in the region of injection which disappears in the course of injection.

⚠ To minimize this risk, the product must be injected very slowly❸.

Once the product has been injected, the needle is left 2 minutes before being withdrawn. **Others** The act must be performed by a medical professional of health familiar with percutaneous punctures of herniated discs.

PRECAUTIONS TO OBSERVE

The viscosity of DISCOGEL® depends on the temperature. Avoid an administration of the product warmed up above room temperature, because gel becomes more liquid and is below optimum viscosity.

To increase its viscosity, DISCOGEL® can be refrigerated just prior to injection.

Administration to children: as the security of this product has not been established for children, extreme caution must be used especially in judging dose.

Do not use the vial or accessories if any of the packaging is damaged since sterility may not be guaranteed.

⚠ ⚠ CAUTION

All the essential elements of kit are in single use. The risk of reusing one of the elements can notably lead to septic risks, evaporation of ethyl alcohol, problems of "syringe-ability" and a defect of effectiveness of the device.

The disc to be treated must represent at least the half of its normal initial height.

CONTRAINDICATIONS

DISCOGEL® is not indicated for patients known to be allergic to one of the components, patients in severe depression or any other condition making the interpretation of pain difficult. DISCOGEL® is not indicated for pregnant woman.

REMARK

The herniated disc is most often linked to a defect of muscle tone and posture. So, it is necessary to follow the treatment with postural re-education. A correctly designed program should avoid both relapse and hernia of other spinal discs.

RECOMMENDATIONS

Simultaneous general or local complementary treatments can be used in case of facet joint pain. Additional, complementary intradiscal therapeutics are also possible in case of treatment of hernias with narrow canal, foraminal, extremely painful and "sleepless" hernias. Depending on the clinical evaluation of the patient, the practitioner may prescribe some painkillers to minimize inflammation.

SIDE EFFECTS

- A theoretical possibility of escape into the periradicular or peridural zone during injection remains possible. This could lead to a feeling of transitional burn in the contact of a root.
- The violent and important high loss of the disc treated, well that unlikely, can lead to severe lumbagoes.
- Transitory radicular irritation with neuropathic pain may appear immediately after the procedure.
- A diffusion of local anesthetic along the introducer needle may induce a transitory sensory-motor deficit.
- Transient pain exacerbation starting one week after DISCOGEL®'s injection may occur.
- Systemic allergic reaction to the substance.
- Discitis and neurologic injury.
- Transient urinary incontinence.
- Sore throat, neck stiffness, dysphagia.

-  **INCOMPATIBILITÄT**
DISCOGEL® should not be blended with other solutions.

CONSERVATION

Keep in a temperature between +2°C and +25°C, protected from light. If this temperature cannot be assured, keep in a refrigerator at 4°C.

COLLECTION AND ELIMINATION OF USED CONSUMABLES

Consumables after use must be discarded according to the rules of good practice.

POST-OPERATIVE ADVICE TO DISPENSE TO THE PATIENTS

The practitioner may advice his/her patient about the post-operative treatment and about the behavior to follow up during the recovery.

	FABRICATEUR GELSCOM S.A.S.	
	8, Avenue Dubna – ZAC CITIS 14200 HÉROUVILLE-SAINT-CLAIR FRANCE	
	contact@gelscom.com www.gelscom.com www.discogel.com	

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